

Point-of-Care Collection Solutions

Simple • Standardized • Small Volumes

Specimen collection requirements for POC are rapidly changing.

Do you have the right tools?

The SARSTEDT S-Monovette® POC Collect Kit, 1 ml Blood Gas Monovette®, and Minivette® POCT offer small sample volumes that meet the requirements for some POC testing cartridges like the i-STAT CHEM8+ cartridge that requires a 95 µl venous sample size.



S-Monovette® POC Collect Kit

- 1.2 ml Lithium-Heparin collection and dispensing tube
- Reduces POC test rejection
- Improved turn-around-time
- Better PBM (Patient Blood Management)
- Closed system = maximum safety
- Compliant with waived testing guidelines

[Brochure Link](#)

[Watch Video](#)



1 ml Blood Gas Monovette®

- Dosed with liquid Ca²⁺ balanced heparin; 50 I.E. heparin/ml of blood
- Suitable for POC blood gas testing
- Can be used with venous or arterial blood
- Better PBM (Patient Blood Management)
- Anaerobic collection with a ventilator cap

[Brochure Link](#)



Minivette® POCT

- Plastic collection and dispensing device
- Drip-free transfer to POC device
- Several volume & additive options from 10 to 200 µl
- Better PBM (Patient Blood Management)
- Closed system = maximum safety
- Compliant with waived testing guidelines

[Brochure Link](#)

[Watch Video](#)

Contact our Point-of-Care Specialist **Cas Wise** for more information!

Email: Lawrence.Wise@Sarstedt.com

Mobile: (828) 461-6530